


Product Datasheet

ERp57 (PDIA3) antibody (orb1334702)

Description	Mouse monoclonal antibody to ERp57 (PDIA3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 140-505 of human PDIA3(NP_005304) produced in E.coli.
Target	PDIA3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	54.2 kDa
Uniprot ID	P30101
NCBI	NM_005313
Dilution Range	WB 1:500, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...