

## Product Datasheet

### C20orf30 (TMEM230) antibody (orb1334738)

**Description** Mouse monoclonal antibody to C20orf30 (TMEM230)

**Species/Host** Mouse

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** FC, WB

**Immunogen** Full length human recombinant protein of human C20orf30(NP\_054864) produced in HEK293T cell.

**Target** TMEM230

**Note** For research use only

**Application notes** For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

**Isotype** IgG2a

**Clonality** Monoclonal


**MW** 13 kDa

**Uniprot ID** [Q96A57](#)

**NCBI** [NM\\_014145](#)

**Dilution Range** WB 1:2000, FLOW 1:100

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells w...