


Product Datasheet

PFDN3 (VBP1) antibody (orb1334765)

Description	Mouse monoclonal antibody to PFDN3 (VBP1)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human VBP1(NP_003363) produced in HEK293T cell.
Target	VBP1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	22.4 kDa
Uniprot ID	P61758
NCBI	NM_003372
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...