


Product Datasheet

PADI4 antibody (orb1334770)

Description	Mouse monoclonal antibody to PADI4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 299-588 of human PADI4(NP_036519) produced in E.coli.
Target	PADI4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	73.9 kDa
Uniprot ID	Q9UM07
NCBI	NM_012387
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...