


Product Datasheet

PPIL6 antibody (orb1334816)

Description	Mouse monoclonal antibody to PPIL6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PPIL6(NP_775943) produced in HEK293T cell.
Target	PPIL6
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	35 kDa
Uniprot ID	Q8IXY8
NCBI	NM_173672
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...