


Product Datasheet

CD31 (PECAM1) antibody (orb1334844)

Description	Mouse monoclonal antibody to CD31 (PECAM1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Immunogen	Full length human recombinant protein of human PECAM1(NP_000433) produced in HEK293T cell.
Target	PECAM1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	82.4 kDa
Uniprot ID	P16284
NCBI	NM_000442
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...