


Product Datasheet

EXOSC3 antibody (orb1334850)

Description	Mouse monoclonal antibody to EXOSC3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human EXOSC3(NP_057126) produced in E.coli.
Target	EXOSC3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	29.4 kDa
Uniprot ID	Q9NQT5
NCBI	NM_016042
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...