

Product Datasheet

Selenium Binding Protein 1 (SELENBP1) antibody (orb1334894)

Description Mouse monoclonal antibody to Selenium Binding Protein 1 (SELENBP1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human SELENB1 (NP_003935) produced in HEK293T cell.

Target SELENBP1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 52.2 kDa

Uniprot ID [Q13228](#)

NCBI [NM_003944](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...