

Product Datasheet

UBE2E3 antibody (orb1334939)

Description	Mouse monoclonal antibody to UBE2E3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human UBE2E3(NP_006348) produced in E.coli.
Target	UBE2E3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22.7 kDa
Uniprot ID	Q969T4
NCBI	NM_006357
Dilution Range	WB 1:1000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...