


Product Datasheet

CRYZL1 antibody (orb1334978)

Description	Mouse monoclonal antibody to CRYZL1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human CRYZL1(NP_665857) produced in HEK293T cell.
Target	CRYZL1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	38.5 kDa
Uniprot ID	O95825
NCBI	NM_145858
Dilution Range	WB 1:500~2000, IF 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...