


Product Datasheet

SERPINE2 antibody (orb1334993)

Description	Mouse monoclonal antibody to SERPINE2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SERPINE2(NP_006207) produced in HEK293T cell.
Target	SERPINE2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	41.8 kDa
Uniprot ID	P07093
NCBI	NM_006216
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...