

Product Datasheet

Transferrin Receptor 2 (TFR2) antibody (orb1335002)

Description Mouse monoclonal antibody to Transferrin Receptor 2 (TFR2)

Species/Host Mouse

Reactivity Human, Monkey, Mouse

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human TFR2(NP_003218) produced in HEK293T cell.

Target TFR2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 88.6 kDa

Uniprot ID [Q9UP52](#)

NCBI [NM_003227](#)

Dilution Range WB 1:500~2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...