


Product Datasheet

HSD17B2 antibody (orb1335023)

Description	Mouse monoclonal antibody to HSD17B2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human HSD17B2(NP_002144) produced in HEK293T cell.
Target	HSD17B2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.6 kDa
Uniprot ID	P37059
NCBI	NM_002153
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...