


Product Datasheet

RBM46 antibody (orb1335037)

Description	Mouse monoclonal antibody to RBM46
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human RBM46(NP_659416) produced in HEK293T cell.
Target	RBM46
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	59.8 kDa
Uniprot ID	Q8TBY0
NCBI	NM_144979
Dilution Range	WB 1:1000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...