

Product Datasheet

Carboxypeptidase A (CPA1) antibody (orb1335057)

Description Mouse monoclonal antibody to Carboxypeptidase A (CPA1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of human CPA1(NP_001859) produced in HEK293T cell.

Target CPA1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 45.5 kDa

Uniprot ID [P15085](#)

NCBI [NM_001868](#)

Dilution Range WB 1:500~2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...