

Product Datasheet

Phosphoribosyl pyrophosphate amidotransferase (PPAT) antibody (orb1335092)

Description Mouse monoclonal antibody to Phosphoribosyl pyrophosphate amidotransferase (PPAT)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, WB

Immunogen Human recombinant protein fragment corresponding to amino acids 42-278 of human PPAT(NP_002694) produced in E.coli.

Target PPAT

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal

MW 57.2 kDa

Uniprot ID [Q06203](#)

NCBI [NM_002703](#)

Dilution Range WB 1:500, FLOW 1:100

Expiration Date 12 months from date of receipt.

