

Product Datasheet

Cystatin C (CST3) antibody (orb1335109)

Description Mouse monoclonal antibody to Cystatin C (CST3)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human CST3(NP_000090) produced in HEK293T cell.

Target CST3

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 13.3 kDa

Uniprot ID [P01034](#)


NCBI [NM_000099](#)

Dilution Range WB 1:2000, IHC 1:150


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts (35ug)...



Western blot analysis of extracts (35ug)...



Western blot analysis of extracts (10ug)...