

Product Datasheet

Cystatin SA (CST2) antibody (orb1335190)

Description Mouse monoclonal antibody to Cystatin SA (CST2)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human CST2(NP_001313) produced in HEK293T cell.

Target CST2

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 14.3 kDa

Uniprot ID [P09228](#)

NCBI [NM_001322](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...