


## Product Datasheet

### GADD34 (PPP1R15A) antibody (orb1335206)

|                            |                                                                                                                                                                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to GADD34 (PPP1R15A)                                                                                                                                                                                                       |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                |
| <b>Reactivity</b>          | Canine, Human, Monkey, Mouse, Rat                                                                                                                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | FC, IF, IHC, WB                                                                                                                                                                                                                                      |
| <b>Immunogen</b>           | Full length human recombinant protein of human PPP1R15A(NP_055145) produced in HEK293T cell.                                                                                                                                                         |
| <b>Target</b>              | PPP1R15A                                                                                                                                                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                |
| <b>Application notes</b>   | For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process. |
| <b>Isotype</b>             | IgG2b                                                                                                                                                                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                           |
| <b>MW</b>                  | 73.3 kDa                                                                                                                                                                                                                                             |
| <b>Uniprot ID</b>          | <a href="#">O75807</a>                                                                                                                                                                                                                               |
| <b>NCBI</b>                | <a href="#">NM_014330</a>                                                                                                                                                                                                                            |
| <b>Dilution Range</b>      | WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100                                                                                                                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                      |



Western blot analysis of HEK293T cells w...