

Product Datasheet

Pallidin (BLOC1S6) antibody (orb1335233)

Description Mouse monoclonal antibody to Pallidin (BLOC1S6)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IHC, WB

Immunogen Full length human recombinant protein of human PLDN(NP_036520) produced in HEK293T cell.

Target BLOC1S6

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 19.6 kDa

Uniprot ID [Q9UL45](#)

NCBI [NM_012388](#)

Dilution Range WB 1:500~2000, IHC 1:150, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...