


Product Datasheet

RNF113B antibody (orb1335344)

Description	Mouse monoclonal antibody to RNF113B
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human RNF113B(NP_849192) produced in HEK293T cell.
Target	RNF113B
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	36.1 kDa
Uniprot ID	Q8IZP6
NCBI	NM_178861
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...