


Product Datasheet

HES1 antibody (orb1335368)

Description	Mouse monoclonal antibody to HES1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length huamn recombinant protein of human HES1(NP_005515) produced in HEK293T cell.
Target	HES1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	29.4 kDa
Uniprot ID	Q14469
NCBI	NM_005524
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...