


Product Datasheet

OBFC2B (NABP2) antibody (orb1335424)

Description	Mouse monoclonal antibody to OBFC2B (NABP2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human OBFC2B(NP_076973) produced in HEK293T cell.
Target	NABP2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22.2 kDa
Uniprot ID	Q9BQ15
NCBI	NM_024068
Dilution Range	WB 1:500~2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...