


Product Datasheet

PDSS2 antibody (orb1335472)

Description	Mouse monoclonal antibody to PDSS2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PDSS2(NP_065114) produced in HEK293T cell.
Target	PDSS2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	43.9 kDa
Uniprot ID	Q86YH6
NCBI	NM_020381
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...