


Product Datasheet

Kv beta 1 (KCNAB1) antibody (orb1335473)

Description	Mouse monoclonal antibody to Kv beta 1 (KCNAB1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human KCNAB1(NP_751891) produced in HEK293T cell.
Target	KCNAB1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	44.5 kDa
Uniprot ID	Q14722
NCBI	NM_172159
Dilution Range	WB 1:1000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...