

Product Datasheet

Kv beta 1 (KCNAB1) antibody (orb1335474)

Description Mouse monoclonal antibody to Kv beta 1 (KCNAB1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Full length human recombinant protein of KCNAB1(NP_751891) produced in HEK293T cell.

Target KCNAB1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 44.5 kDa

Uniprot ID [Q14722](#)

NCBI [NM_172159](#)

Dilution Range WB 1:1000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...