


Product Datasheet

Adiponectin (ADIPOQ) antibody (orb1335524)

Description	Mouse monoclonal antibody to Adiponectin (ADIPOQ)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 19-244 of human ADIPOQ(NP_004788) produced in HEK293T cell.
Target	ADIPOQ
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	24.9 kDa
Uniprot ID	Q15848
NCBI	NM_004797
Dilution Range	IF 1:100, FLOW 1:100, WB 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...