


Product Datasheet

katanin p80 (KATNB1) antibody (orb1335597)

Description	Mouse monoclonal antibody to katanin p80 (KATNB1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human KATNB1(NP_005877) produced in HEK293T cell.
Target	KATNB1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	72.2 kDa
Uniprot ID	Q9BVA0
NCBI	NM_005886
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...