


Product Datasheet

Cytochrome P450 2C9 (CYP2C9) antibody (orb1335608)

Description	Mouse monoclonal antibody to Cytochrome P450 2C9 (CYP2C9)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human CYP2C9(NP_000762) produced in HEK293T cell.
Target	CYP2C9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	55.4 kDa
Uniprot ID	P11712
NCBI	NM_000771
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...