


Product Datasheet

APBB3 antibody (orb1335639)

Description	Mouse monoclonal antibody to APBB3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human APBB3(NP_573419) produced in HEK293T cell.
Target	APBB3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	52.4 kDa
Uniprot ID	O95704
NCBI	NM_133173
Dilution Range	WB 1:500, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...