


Product Datasheet

TBC1D21 antibody (orb1335647)

Description	Mouse monoclonal antibody to TBC1D21
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length recombinant protein of human TBC1D21(NP_699187) produced in HEK293T cell.
Target	TBC1D21
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	39 kDa
Uniprot ID	Q8IYX1
NCBI	NM_153356
Dilution Range	WB 1:500~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...