


Product Datasheet

ALS2CR1 (NIF3L1) antibody (orb1335668)

Description	Mouse monoclonal antibody to ALS2CR1 (NIF3L1)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human NIF3L1(NP_068596) produced in HEK293T cell.
Target	NIF3L1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	38.8 kDa
Uniprot ID	Q9GZT8
NCBI	NM_021824
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...