


Product Datasheet

PIG3 (TP53I3) antibody (orb1335705)

Description	Mouse monoclonal antibody to PIG3 (TP53I3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TP53I3(NP_671713) produced in HEK293T cell.
Target	TP53I3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.4 kDa
Uniprot ID	Q53FA7
NCBI	NM_147184
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...