


Product Datasheet

Acetyl CoA synthetase (ACSS2) antibody (orb1335737)

Description	Mouse monoclonal antibody to Acetyl CoA synthetase (ACSS2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ACSS2(NP_061147) produced in HEK293T cell.
Target	ACSS2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	78.4 kDa
Uniprot ID	Q9NR19
NCBI	NM_018677
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...