


Product Datasheet

ELK3 antibody (orb1335744)

Description	Mouse monoclonal antibody to ELK3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ELK3(NP_005221) produced in HEK293T cell.
Target	ELK3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	44.1 kDa
Uniprot ID	P41970
NCBI	NM_005230
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...