


Product Datasheet

PPAR alpha (PPARA) antibody (orb1335797)

Description	Mouse monoclonal antibody to PPAR alpha (PPARA)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 145-420 of human PPARA(NP_005027) produced in E.coli.
Target	PPARA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	52 kDa
Uniprot ID	Q07869
NCBI	NM_005036
Dilution Range	WB 1:200~500, IF 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent
staining of COS7
cells...