


Product Datasheet

FDFT1 antibody (orb1335842)

Description	Mouse monoclonal antibody to FDFT1
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human FDFT1(NP_004453) produced in HEK293T cell.
Target	FDFT1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	47.9 kDa
Uniprot ID	P37268
NCBI	NM_004462
Dilution Range	WB 1:200~500, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...