


Product Datasheet

LGALS3BP antibody (orb1335851)

Description	Mouse monoclonal antibody to LGALS3BP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment cprresponding to amino acids 19-300 of human LGALS3BP(NP_005558) produced in E.coil.
Target	LGALS3BP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	65.2 kDa
Uniprot ID	Q08380
NCBI	NM_005567
Dilution Range	WB 1:2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...