


Product Datasheet

TMPRSS5 antibody (orb1335887)

Description	Mouse monoclonal antibody to TMPRSS5
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human TMPRSS5(NP_110397) produced in HEK293T cell.
Target	TMPRSS5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	49.4 kDa
Uniprot ID	Q9H3S3
NCBI	NM_030770
Dilution Range	WB 1:500, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...