


Product Datasheet

CARKL (SHPK) antibody (orb1335948)

Description	Mouse monoclonal antibody to CARKL (SHPK)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human SHPK(NP_037408) produced in HEK293 cell.
Target	SHPK
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	51.3 kDa
Uniprot ID	Q9UJH6
NCBI	NM_013276
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...