

Product Datasheet

TRAP alpha (SSR1) antibody (orb1335992)

Description Mouse monoclonal antibody to TRAP alpha (SSR1)

Species/Host Mouse

Reactivity Canine, Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human SSR1(NP_003135) produced in HEK293T cell.

Target SSR1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 32.1 kDa

Uniprot ID [P43307](#)

NCBI [NM_003144](#)

Dilution Range WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...