

Product Datasheet

XPF (ERCC4) antibody (orb1335994)

Description	Mouse monoclonal antibody to XPF (ERCC4)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human ERCC4(NP_005227) produced in HEK293 cell.
Target	ERCC4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	104.3 kDa
Uniprot ID	Q92889
NCBI	NM_005236
Dilution Range	WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.

