


Product Datasheet

DHFRL1 (DHFR2) antibody (orb1335995)

Description	Mouse monoclonal antibody to DHFRL1 (DHFR2)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 4-185 of human DHFRL1(NP_789785) produced in E.coli.
Target	DHFR2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	21.4 kDa
Uniprot ID	Q86XF0
NCBI	NM_176815
Dilution Range	WB 1:500~2000, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...