


Product Datasheet

Tropomodulin 1 (TMOD1) antibody (orb1336080)

Description	Mouse monoclonal antibody to Tropomodulin 1 (TMOD1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human TMOD1(NP_003266) produced in HEK293 cell.
Target	TMOD1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	40.4 kDa
Uniprot ID	P28289
NCBI	NM_003275
Dilution Range	WB 1:500~2000, IHC 1:50, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...