

Product Datasheet

SM22 alpha (TAGLN) antibody (orb1336127)

Description Mouse monoclonal antibody to SM22 alpha (TAGLN)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full length human recombinant protein of human TAGLN(NP_003177) produced in HEK293T cell.

Target TAGLN

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 22.4 kDa

Uniprot ID [Q01995](#)

NCBI [NM_003186](#)

Dilution Range WB 1:500~2000, IHC 1:150, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...