


Product Datasheet

MDS028 (ITFG2) antibody (orb1336162)

Description	Mouse monoclonal antibody to MDS028 (ITFG2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human ITFG2 (NP_060933) produced in HEK293T cell.
Target	ITFG2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	49.1 kDa
Uniprot ID	Q969R8
NCBI	NM_018463
Dilution Range	WB 1:500, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...