

Product Datasheet

LIM Kinase 1 (LIMK1) antibody (orb1336193)

Description Mouse monoclonal antibody to LIM Kinase 1 (LIMK1)

Species/Host Mouse

Reactivity Canine, Human, Monkey, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full length human recombinant protein of human LIMK1(NP_002305) produced in HEK293T cell.

Target LIMK1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 72.4 kDa

Uniprot ID [P53667](#)

NCBI [NM_002314](#)

Dilution Range WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...