


Product Datasheet

ALDH1A3 antibody (orb1336330)

Description	Mouse monoclonal antibody to ALDH1A3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-100 and 413-512 of human ALDH1A3 (NP_000648) produced in E.coli.
Target	ALDH1A3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	55.9 kDa
Uniprot ID	P47895
NCBI	NM_000693
Dilution Range	WB 1:200~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...