


Product Datasheet

MAB21L3 antibody (orb1336423)

Description	Mouse monoclonal antibody to MAB21L3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human C1orf161 (NP_689580) produced in HEK293T cell.
Target	MAB21L3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	42.2 kDa
Uniprot ID	Q8N8X9
NCBI	NM_152367
Dilution Range	WB 1:2000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...