


Product Datasheet

MEF2C antibody (orb1336424)

Description	Mouse monoclonal antibody to MEF2C
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human MEF2C(NP_002388) produced in HEK293T cell.
Target	MEF2C
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	51 kDa
Uniprot ID	Q06413
NCBI	NM_002397
Dilution Range	WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...